

# GUS® DISINFECTION SOAK STATIONS FOR SMALL SCOPES WALL MOUNTED SERIES

## FEATURES:

- Safely eliminates fumes from Glutaraldehyde and OPA
- Probes are suspended in containers preventing damage to the delicate tips
- Fumes are drawn away from the operator and into a patented carbon filter where they are neutralized
- Locking safety door allows you to install the system in a patient area
- Meets OSHA and JCAHO requirements

## SPECIFICATIONS:

Model	G10ENT	G14ENT-3
Number of scopes soaked	1	2
Face Velocity (FPM)	45	45
Electrical (115V 50/60Hz)	0.5A	0.5A
<b>Exterior Dimensions</b>		
Width	10.5"	15.5"
Depth	10.5"	10.5"
Height	32" (with tubes)	32" (with tubes)
Weight (lbs)	24	28

\*Allow 2" additional space on right for the power cord and 1.5" above the unit for the clean air exhaust.  
UL and CE approved. 115/230 V models available.



**GUS Model G10ENT:**  
Two 19" tubes, one for soak  
and one for rinse

## OPTIONS:



**Glute Out® Neutralizer**



**Spill Kits**



**Trays**



**Scope Holders**



**GUS Model G14ENT-3:**  
Three 19" tubes, two for soak  
and one for rinse



PCI Medical | 6 Winter Ave | Deep River, CT 06417

Tel: USA (800) 862-3394 | International (860) 526-2862 | Fax: (860) 526-3081 | info@pcimedical.com