

# GUS® DISINFECTION SOAK STATIONS FOR LARGE SCOPES SPACE SAVER SERIES

## FEATURES:

- Meets Joint Commissions & OSHA requirements
- Two soak/rinse trays in various sizes
- Removes vapors & protects staff
- No ductwork or plumbing required
- Built-in chemical transfer pump
- Built-in GFCI

The GUS® Space Saver Series fit under an upper wall cabinet. They have the same features as the Elevated Series, but are only 16" high.

Three widths available: 30", 36" and 40"



*GUS Model G32-S  
Shown Above*

## SPECIFICATIONS:

Model	G25-S	G32-S	G36-S
Face Velocity (FPM) (lower door open & trays in place)	58	50	70
Electrical (115V 50/60Hz)	1.0A	1.0A	2.0A
<b>Exterior Dimensions*</b>			
Width	30"	36"	40"
Depth	27"	27"	27"
Height	16.5"	19"	19"
Weight (lbs)	85	105	140

\*Allow 2" additional space on right for the power cord and 1.5" above the unit for the clean air exhaust.

UL and CE approved. 115/230 V models available.

## OPTIONS:



*Glute Out® Neutralizer*



*Spill Kits*



*Trays*



*Optional Tray Inserts*



PCI Medical | 6 Winter Ave | Deep River, CT 06417

Tel: USA (800) 862-3394 | International (860) 526-2862 | Fax: (860) 526-3081 | info@pcimedical.com