

# GUS® DISINFECTION SOAK STATIONS FOR LARGE SCOPES ELEVATED SERIES

## FEATURES:

- Meets Joint Commissions & OSHA requirements
- Two soak/rinse trays in various sizes
- Removes vapors & protects staff
- No ductwork or plumbing required
- Built-in chemical transfer pump
- Built-in GFCI

The GUS® Elevated Series gives extra working height for an easier disinfection process.

An automatic two-speed blower increases airflow when the lower door is open. When the door is closed the unit is in its energy saving mode. The built-in disposal system allows for safe and hands free transfer of used disinfectant. Three widths available: 27", 33.25" and 40"



**GUS Model G32-E  
Shown Above**

## SPECIFICATIONS:

Model	G25-E	G32-E	G36-E
Face Velocity (FPM) (lower door open & trays in place)	58	50	70
Electrical (115V 50/60Hz)	1.0A	1.0A	2.0A
<b>Exterior Dimensions*</b>			
Width	27"	33.25"	40"
Depth	26.5"	26.5"	26.5"
Height	25"	25"	25"
Weight (lbs)	90	110	150

\*Allow 2" additional space on right for the power cord and 1.5" above the unit for the clean air exhaust.

UL and CE approved. 115/230 V models available.

## OPTIONS:



**Glute Out® Neutralizer**



**Spill Kits**



**Trays**



**Optional Carts**



PCI Medical | 6 Winter Ave | Deep River, CT 06417

Tel: USA (800) 862-3394 | International (860) 526-2862 | Fax: (860) 526-3081 | info@pcimedical.com