

GUS® DISINFECTION SOAK STATIONS FOR SMALL SCOPES COUNTER-TOP SERIES

FEATURES:

- Meets OSHA & JCAHO requirements
- Protects scopes during soak & rinse
- Allows for safe transfer of used disinfectant
- Removes vapors & protects staff
- No ductwork or plumbing required
- Instrument Specific

SPECIFICATIONS:

Model	G17	G18-E
Face Velocity (FPM) (lower door open & trays in place)	45	70
Electrical (115V 50/60Hz)	1.0A	1.0A
Exterior Dimensions		
Width	18.25"	19.00"
Depth	24.5"	27.5"
Height	19.5"	24.5"
Weight (lbs)	60	70

*Allow 2" additional space on right for the power cord and 1.5" above the unit for the clean air exhaust.
UL and CE approved. 115/230 V models available.



The G17 is designed to fit under an upper wall cabinet.
* Comes with a stand-alone disposal pump (DS02)
* Manual high/low air speed switch.

OPTIONS:



Glute Out® Neutralizer



Spill Kits



Trays



Optional Tray Inserts



The G18-E gives optimal working height space.
* Built-in disposal system that allows for safe transfer of used disinfectant
* Automatic two-speed blower switches to high when door is opened, increasing airflow through patented carbon filter.



PCI Medical | 6 Winter Ave | Deep River, CT 06417

Tel: USA (800) 862-3394 | International (860) 526-2862 | Fax: (860) 526-3081 | info@pcimedical.com